

$TamiSolve^{\rm TM}$ NxG as a suitable solvent for cost-efficient green SPPS

Jan Pawlas,1 Mikael Nilsson, Biljana Antonic, Thomas Svensson, Jens Finnman, Jon H. Rasmussen

PolyPeptide Group, Limhamnsvägen 108, PO BOX 30089, 20061 Limhamn, Sweden

In recent years replacing hazardous solvents by greener alternatives has become a point of emphasis in chemical and pharma industries.² For example, in peptide synthesis great advances have been made in the past decades although it has been stated that more attention needs to be paid to the impact of solvents and reagents, used in this field, on health and environment.3 Specifically, solid-phase peptide synthesis (SPPS) is a dominant synthetic methodology which has traditionally relied on use of large amounts of N,N-dimethyl formamide (DMF), a polar aprotic solvent with reprotoxic⁴ and potentially carcinogenic characteristics.⁵ Consequently, SPPS in which DMF is substituted for less harmful alternatives has become an active area of research in recent years.⁶ Our aim with green SPPS is to put forth protocols suitable for large scale synthesis of pharmaceutical peptides that would not only replace DMF with green solvents, but improve the synthesis at the same time. To this end we set out to evaluate Nbutylpyrrolidinone (NBP, cas nr. 3470-98-2) also known as TamiSolveTM NxG⁷ as a solvent for green SPPS. NBP has been reported as a biodegradable alternative for DMF in several common synthetic transformations⁸ and our interest in NBP was triggered by the fact that in NBP, Fmoc-Cys(Trt)-OH in the presence of coupling additive Oxyma was completely stable. This is a remarkable improvement over DMF, in which Fmoc-Cys(Trt)-OH/ Oxyma decomposes extensively, conceivably via an Oxyma N-oxyl radical induced mechanism (Figure 1).9 During the course of our NBP evaluation a report on replacing DMF with NBP (1:1) in SPPS of Octreotide has been disclosed, in which NBP gave slightly inferior purity of the peptide (80 %) compared to DMF (86%).10 Unfortunately, NBP as a stand alone SPPS solvent is too expensive to be considered suitable for cost-efficient peptide synthesis. We



therefore set out to examine NBP in SPPS together with less polar, inexpensive green cosolvent(s). We envisioned that this approach would result in inexpensive solvent systems well suited for SPPS on PS resins commonly used in large scale manufacturing. As a model substrate we chose the Aib-ACP 10 mer 1.

Figure 1. HPLC of 0.1M Fmoc-Cys(Trt)-OH/Oxyma in DMF (in black) and in NBP (in blue) after 72 h at rt. For identities of Cys degradants in the DMF sample see ref 6a.

which in a recent green synthesis on PS resin gave the target molecule in 54% purity.⁶⁰

H-Val¹-Gln²-Aib³-Aib⁴-Ile⁵-Asp⁶-Tyr⁷-Ile⁸-Asn⁹-Gly¹⁰-NH₂(1)

We chose to examine the most difficult part of the Aib-ACP SPPS (Val1-Aib4) for which we used a Fmoc-5-10 resin prepared by our earlier green SPPS protocol relying on the use of DMSO/EtOAc as the solvent system.5a Throughout this study we used 4-methylpiperidine (4-MP) for all de'Fmocings, Oxyma/DIC for all couplings and mild conventional heating (40 - 45 °C) to attain suitable SPPS kinetics.11

In our first experiment we examined four coupling solvents i) DMF ii) NBP/EtOAc (1:9) and iv) DMSO/EtOAc (1:9)

Table 1. SPPS of Aib-ACP using different coupling solvents			
Fmoc-5-1	0-RMG-AMS	MP in DMF (2 x 15 min) mL DMF wash (5 min ea uuv AA/Oxyma/DIC (1:1 : mL DMF wash (5 min)	
150 mg, 0	5. TFA/T 6. ether	IS/H ₂ O (92.5:5.0:2.5), 1 h precipitation	
Entry	Coupling solvent	LC-HRMS purity (?	
1	DMF	75.2	
0			
2	DMSO/EtOAc, 1:9	84.9	
3	DMSO/EtOAc, 1:9	84.9 76.4	
	Table 1. 5 Fmoc-5-1 150 mg, 0 Entry 1	Table 1. SPPS of Aib-ACP using diff 1. SPA of Aib-ACP using diff	

Emoc-5-10-RMG- 50 mg, 0.21 M, 0.032 mmol 50 mg, 0.21 M, 0.032 mmol 51 M = 1 M =				
Entry	Coupling solvent	LC-HRMS purity (%)	des Aib (%)	
1	DMF	75.2	12.8	
2	DMSO/EtOAc, 1:9	84.9	2.6	
3	NBP	76.4	10.4	
4	NBP/EtOAc, 1:9	83.6	1.4	

SPPS than DMF (entry 1) and NBP (entry 3) do. We next examined solvents for Fmoc deprotections, using NBP/EtOAc as the solvent for couplings but changing the NBP/EtOAc ratio from 1:9 to 1:1 to ensure suitable solubility of all common AA derivatives at sufficient concentrations (>~0.2 M). We tested DMF, NBP, NRP/EtOAc

And Control and Con				 (1:1) and NBP/EtOAc (1:4) as de'Fmocing solvents (Table 	
Entry	de'Fmocing solvent	LC-HRMS purity (%)	des Aib (%)	2) and found	
1	DMF	77.7	4.8	that both	
2	NBP	81.0	5.3	NBP/EtOAc	
3	NBP/EtOAc, 1:1	83.0	2.1	mixtures	
4	NBP/EtOAc, 1:4	82.5	2.5	(entries 3 and 4) gave higher	

purities of Aib-ACP than DMF (entry 1) and NBP (entry 2) did. It is worth noting that the crude peptides from the

EtOAc. In an attempt to suppress the Ac-4-10-NH2 impurity we reexamined solvents for couplings, in which we tested 3 green 4-10-NH, DMF in de'Fmocings Aib-ACP des Aib coupling solvent mixtures NBP/EtOAc. NBP/i-PrOAc NBP in de'Fmocings and NBP/MeCN respectively (Table 3). NBP/EtOAc (1:1) in de'Fmocings Importantly, AcOH content in all solvents in NBP/EtOAc (1:4) in de'Fmocings this experiment was determined

NBP/EtOAc runs did not contain any peaks which were not present in the crudes using DMF and NBP as solvents

(Figure 2). On the other hand, all four crudes contained comparable amounts of Ac-4-10-NH₂ (3 - 4 %),

conceivably due to the background acetylation during the Aib³ coupling caused by the presence of AcOH in



ppm. We used NBP/EtOAc 1:1 as de'Fmocing solvent throughout and found that NBP/i-PrOAc and NBP/MeCN as coupling solvents gave comparable purity of Aib-ACP (entries 2 - 3),12 both of which were higher than the purity obtained with NBP/EtOAc (entry 1). Compared to the amounts of Ac-4-10-NH₂ seen in Figure 2 the content of this Ac truncation in the three crudes in Table 3 was significantly reduced, down to the levels which can be explained by the presence of minute amounts of AcOH in AA derivatives.¹³ Finally, having found that NBP/EtOAc constitutes a suitable

solvent for Table 3. SPPS of Aib-ACP using different coupling solvents II de'Fmocings and NBP/i-PrOAc or NBP/MeCN are suitable solvents for couplings we set out to examine green solvent wash protocols. As washes

Fmoc.5-10-RMG-AM 4.8% 4-MP in NBP/EIOAc,1:1 (2 x 15 min) 2.1 x 10 mL DMF wash (5 min) 2.0 equiv AXOyma/DIC (1:1:2) in solvent (25 min) 4.1 x 10 mL DMF wash (5 min) 4.5 x 10 mL D				
Entry	Coupling solvent	LC-HRMS purity (%)	des Aib (%)	
1	NBP/EtOAc, 1:1	83.0	2.5	
2	NBP/iPrOAc, 1:1	88.8	1.9	
3	NBP/MeCN, 1:1	88.5	1.9	

account for the vast majority of solvent consumption in SPPS we chose not to investigate NBP, which would be too costly if used in the large amounts needed for the washes. Instead, we attempted to i) omit the washes after couplings ii) use EtOAc instead of DMF in washes after de'Fmocings (Table 4) to render the washing as simple and inexpensive as possible.¹⁴ These experiments showed that solvent washes after couplings can be fully eliminated without any erosion in purity of Aib-ACP (entries 1 and 2 vs 3 and 4). On the other hand, replacing DMF with EtOAc as solvent after de'Fmocings resulted in a slight decrease in Aib-ACP purity (entries 1 and 3 vs 2 and 4). In fact, this observation is in a good agreement with our previous evaluation of green solvent washes after de'Fmocings in which we determined that neat EtOAc is inferior to DMF whereas using 2% DMSO in EtOAc orobl

Table 4. SPPS of Aib-ACP using different solvent wash protocols $\begin{pmatrix} 1.8% - 4MP \text{ in MBP/EIOAc.1.1 (2 x 15 min)} \\ 2. 1 x 10 mt. solvent wash (5 min) \\ 3. 2.0 equiv AAOxyma/DIC (11:3) in NBP/MeCN, 1:1 (25 min) \\ 4. 1 x 10 mt. solvent wash (5 min) \end{cases}$					e crude purities as with DMF
Finded for Render FFATTIS/H₂O (92.5:5.0:2:5), 1 h at rt H-1-10-NH₂ 150 mg, 0.21 M, 0.032 mmol 6. ether precipitation (Aib-ACP)				can be achieved. ⁶	
Entry Sol after	vent wash de'Fmocing	Solvent wash after coupling	LC-HRMS purity (%) ²	des Aib (%)	It is worth
1	DMF	DMF	89.5	3.2	stating that
2	EtOAc	DMF	84.4	2.5	86.3% attained in
3	DMF	none	90.5	1.2	Table 4,
4	EtOAc	none	86.3	1.9	entry 3

synthesis is to our knowledge the highest purity of Aib-ACP synthesized by a fully green SPPS protocol on a PS resin

In summary we examined TamiSolveTM NxG (NBP) as a solvent for green SPPS and found that NBP with inexpensive, common green solvents such as EtOAc and MeCN constitutes a highly competent solvent system which in Aib-ACP SPPS performs better than both DMF and NBP. As a solvent for washes after chemical steps in SPPS we propose to omit washes after couplings and use EtOAc with small amount of a polar cosolvent such as DMSO or NBP for washes after de'Fmocings.^{6a} This green SPPS approach combines the use of solvents which work well in the chemical steps of the synthesis with inexpensive solvent washing, rendering the methodology well suited for cost-efficient peptide manufacturing

to be quite negligible, below the LOO of 50

Correspondence to jap@polypeptide.com. See e.g. Byrne, F. P., Jin, S.: Paggiola, G.: Petchey, T. H. M.; Clark, J. H.; Farmer, T. J.; Hunt, A. J.; McElroy, C. R.; Sherwood, J. Sustain. Chem. Process. 2016, 4, 1. Data, S.; Sood, A.: Jrönk, M. Current Org, Synth. 2011, 8, 262. https://www.echa.europa.eudocuments/1016/23/f309004e393-4th3-8aa-6664e7999e0e66 https://moorphic.inst.frESG/Pholicianos/interrupti-10402.pdf See for example a) Pawlas, J.; Lundqvis, M.; Svensson, T.; Nilsson, M.; Rasmusen, J. H.; Poster P-280 at 25th American Peptide Symposium, June 2017, Whisler, B See for example a) Pawlas, J.; Lundqvis, M.; Svensson, T.; Nilsson, M.; Rasmusen, J. H.; Poster P-280 at 25th American Peptide Symposium, June 2017, Whisler, B ¹and/monographs intre (IESQFDablications/intermerp/14-002.pdf
¹and/monographs inter (IESQFDablications/intermerp/14-002.pdf)
¹And/monographs inter (IESQFDablications/intermerp/14-002.